

=> d his

(FILE 'HOME' ENTERED AT 17:53:50 ON 04 JUN 2008)

FILE 'CA' ENTERED AT 17:54:33 ON 04 JUN 2008

L1 6887 S LUMINOL OR ISOLUMINOL

L2 286427 S (BLOOD OR HEMOGLOBIN OR HEMIN OR PEROXIDASE) (5A) (DETECT? OR DETERMIN? OR MEASUR? OR MONITOR? OR FIND? OR SEARCH? OR ASSAY? OR ANALY? OR SENSE# OR SENSING)

L3 1147 S L1 AND L2

L4 2570 S (HYDROGEN PEROXIDE OR H2O2 OR HOOH) AND L1

L5 369 S (NAOH OR KOH OR (NA OR SODIUM OR POTASSIUM OR K) (1A) (OH OR HYDROXIDE)) AND L1

L6 170 S L4 AND L5

L7 166 S L6 AND (CHEMILUM? OR LUMINESC? OR EMISSION OR EMITTED OR EMITTING)

L8 21 S L2 AND L6

L9 146 S L7 NOT L8

L10 11 S L9 AND (MM OR MMOL?)

L11 32 S L8, L10

=> d l11 1-32 bib, ab, it

L11 ANSWER 6 OF 32 CA COPYRIGHT 2008 ACS on STN

AN 139:347002 CA

TI Composition and process for detecting and locating traces of blood

IN Lefebvre Despeaux, Jean Marc; Blum, Loic

PA Roc Import, Monaco

SO Fr. Demande, 29 pp.

PI FR 2839155 A1 20031031 FR 2002-5230 20020425

US 20050176082 A1 20050811 US 2004-509326 20041006

PRAI FR 2002-5230 A 20020425

AB The title compn. comprises luminol or its deriv., hydrogen peroxide and NaOH, dild. in a solvent, preferably water. The compn. is usable in forensics or for detecting traces of game blood during hunt.

L11 ANSWER 9 OF 32 CA COPYRIGHT 2008 ACS on STN

AN 137:72183 CA

TI On-line determination of transition metal ions by low-pressure ion chromatography with chemiluminescence detection

AU Zhou, Guang-Ming; Liu, Ling; Yang, Guang-Ming; Zhang, Xin-Shen

CS School of Chemistry and Chemical Engineering, Southwest Normal University, Chungking, 400715, Peop. Rep. China

SO Sepu (2002), 20(3), 265-268

LA Chinese

AB The sepn. and detn. of five transitional metal ions, Co<sup>2+</sup>, Cu<sup>2+</sup>, Mn<sup>2+</sup>, Fe<sup>2+</sup> and Cr<sup>3+</sup> by low-pressure ion chromatog. with chemiluminescence detection were investigated for the first time. The sepn. was carried out on a C3 ion chromatog. column by using a step gradient elution with the mobile phase of oxalic acid-citric acid soln. The five transitional-metal ions could be sepd. simultaneously. Luminol-H<sub>2</sub>O<sub>2</sub>-Mn<sup>+</sup> system was used as the flow injection chemiluminescence anal. system. Under the optimum conditions with 0.5 mmol L<sup>-1</sup> luminol, 0.05 mol L<sup>-1</sup> H<sub>2</sub>O<sub>2</sub> and 0.020 mol L<sup>-1</sup> NaOH, the linear ranges (mg L<sup>-1</sup>) for Co<sup>2+</sup>, Cu<sup>2+</sup>, Mn<sup>2+</sup>, Fe<sup>2+</sup> and Cr<sup>3+</sup> were 0.001 - 0.1, 0.1 - 6, 0.06 - 4, 0.03 - 5 and

0.025 - 1 resp.; the detection limit ( $\mu\text{g L}^{-1}$ ) for  $\text{Co}^{2+}$ ,  $\text{Cu}^{2+}$ ,  $\text{Mn}^{2+}$ ,  $\text{Fe}^{2+}$  and  $\text{Cr}^{3+}$  were 0.85, 85, 42, 21 and 20. The present method is simple, rapid, accurate and has been successfully applied to detect  $\text{Co}^{2+}$  in vitamin B12 and  $\text{Co}^{2+}$ ,  $\text{Cu}^{2+}$ ,  $\text{Mn}^{2+}$ ,  $\text{Fe}^{2+}$  and  $\text{Cr}^{3+}$  in human hair sample with satisfactory results.

L11 ANSWER 17 OF 32 CA COPYRIGHT 2008 ACS on STN

AN 119:43857 CA

OREF 119:7863a,7866a

TI Chemiluminescence determination of horseradish peroxidase and its conjugates with several systems

AU Yang, Xiucen; Wu, Liping; Deng, Anping; Chen, Dajing

CS Res. Lab. Appl. Chem., West China Univ. Med. Sci., Chengdu, 610041, Peop. Rep. China

SO Analytical Letters (1993), 26(6), 1065-71

AB Four systems used for chemiluminescence detn. of horseradish peroxidase [luminol- $\text{H}_2\text{O}_2$  (NaOH), pyrogallol- $\text{H}_2\text{O}_2$  (pH 6.5), luminol- $\text{H}_2\text{O}_2$ -para-iodophenol (pH 8.5), and luminol- $\text{H}_2\text{O}_2$ -para-hydroxybiphenyl (pH 8.5)] were studied. The exptl. conditions were optimized and linear dynamic ranges and limits of detection detd. The HRP conjugates with different antibodies were also detd., on the basis of which methods for detn. of insulin, gentamycin, and Salmonella anatis by chemiluminescent immunoassay were established.

=> log y

STN INTERNATIONAL LOGOFF AT 18:08:43 ON 04 JUN 2008